1 Comparison of the elastic recovery and strain-in-

- 2 compression of commercial and novel vinyl polysiloxane
- 3 impression materials incorporating a novel crosslinking
- 4 agent and a surfactant
- 5
- 6 Shahab Ud Din ¹, Osama Khattak ^{2*}, Farooq Ahmad Chaudhary ^{1*}, Asfia Saeed ³, Azhar Iqbal ²,
- 7 Jamaluddin Syed ⁴, Alaa Ahmed Kensara ⁵, Thani Alsharari ⁶, Mohammed Mustafa ⁷, Sherif
- 8 Elsayed Sultan ⁸, Mangala Patel ⁹
- 9 ¹ School of Dentistry (SOD), Shaheed Zulfiqar Ali Bhutto Medical University (SZABMU),
- 10 Islamabad, Pakistan.
- 11 ² Department of Restorative Dentistry, Jouf University, Sakaka, Saudi Arabia.
- ³ Shifa College of Dentistry, Shifa Tameer-e-Millat University, Islamabad, Pakistan.
- ⁴ Oral Basic Clinical Sciences Department, faculty of Dentistry, King Abdulaziz university,
 Jeddah, Saudi Arabia.
- ⁵ Department of Oral and Maxillofacial Prosthodontics, King Abdul Aziz University Dental
- 16 Hospital, , Jeddah, Saudi Arabia.
- ⁶ Department of Restorative and dental science, faculty of dentistry, Taif University, Taif,
- 18 11099, Saudi Arabia.
- ⁷ Department of Conservative Dental Sciences, College of Dentistry, Prince Sattam Bin
 Abdulaziz University, Al-Kharj 11942, Saudi Arabia
- ⁸ Department of fixed prosthodontics, Tanta University, Tanta, Egypt.
- ⁹ Centre for Oral Bioengineering (Dental Physical Sciences Unit), Bart's and The London,
- 23 School of Medicine and Dentistry, Queen Mary University of London, UK.
- 24
- 25 Corresponding Authors:
- 26 1. Osama Khattak
- 27 Department of Restorative Dentistry, Jouf University, Sakaka, Saudi Arabia.
- 28 Email address: dr.osama.khattak@jodent.org
- 29 2. Farooq Ahmad Chaudhary

30 School of Dentistry (SOD), Shaheed Zulfiqar Ali Bhutto Medical University (SZABMU),

31 Islamabad, Pakistan.

32 Email: chaudhary4@hotmail.com

33 Abstract

34 This study aims to formulate experimental vinyl polysiloxane (VPS) impression materials and 35 compare their elastic recovery and strain-in-compressions with three commercial VPS materials 36 (Aquasil, Elite, and Extrude). Five experimental materials (Exp), two hydrophobic (Exp-I and II) 37 and three hydrophilic (Exp-III, IV and V) were developed. Exp 1 contained vinyl-terminated polydimethyl siloxane and a conventional cross-linking agent (poly methylhydrosiloxane), while Exp-38 39 II contained a novel cross-linking agent that is tetra-functional dimethyl-silyl-ortho-silicate 40 (TFDMSOS). Exp III-V (hydrophilic materials) were formulated by incorporating different 41 concentrations of non-ionic surfactant (Rhodasurf CET-2) into Exp II formulation. Measurement 42 of elastic recovery and strain-in-compression for commercial and experimental materials were 43 performed according to ISO4823 standard using the calibrated mechanical testing machine (Tinius 44 Olsen). One-way Analysis of Variance (one-way ANOVA) and Tukey's posthoc (HSD) test were 45 used for statistical analysis and a p-value of ≤ 0.05 was considered significant. Exp-I has 46 statistically similar values to commercial VPS. The Exp-II showed the highest elastic recovery, 47 while % elastic recovery was reduced with the addition of the non-ionic surfactant (Rhodasurf 48 CET-2). The % reduction was directly related to the concentration of Rhodasurf CET-2. In 49 addition, Exp II had significantly higher strain-in-compression values compared to Exp-I and 50 commercial materials. These values were further increased with the addition of a non-ionic 51 surfactant (Rhodasurf CET-2) was added (Exp-III, IV and V).

52

53 Introduction

Dental impression refers to a negative imprint of oral hard and soft tissues ¹. An accurate impression is of utmost importance for the successful fabrication of a prosthesis. Ideally, the impression material should have good wettability, accuracy, elasticity and minimal distortion on removal and storage ¹. Impression materials are compressed against the tray, especially while recording undercut areas, and on the removal of impression from the mouth. The degree of distortion of the material depends on the severity of the undercut, elastic recovery of the material, the time the material is kept in the compressed state and storage conditions ^{1,2}.

61

The elastic recovery of impression material is the capacity of the material to revert to its original position, without significant distortion after being strained, when the deforming force is removed (Hamalian et al. 2011). It is due to the presence of folded polymeric segments within the material, which coil and uncoil during loading and unloading. Therefore, the greater the elastic recovery of the material, the more precise the prosthesis.

67

The likelihood of permanent deformation increases on slow removal of an impression as the 68 material is stressed for longer duration (Balkenhol et al. 2010; Din et al. 2018b; Din et al. 2021; 69 70 Din et al. 2022; Hondrum 1994; Mandikos 1998). None of the impression materials has 100% 71 elastic recovery(Hamalian et al. 2011), rather most elastomeric materials exhibit time-dependent 72 recovery from deformation (viscoelasticity) (Braden et al. 1997; Darvell 2009; Goldberg 1974). 73 The elastic recovery of these materials depends on their composition, such as the pre-polymer, 74 cross-linking agents, and fillers (Deb 1998; Din et al. 2018a; Din et al. 2018b; Klooster et al. 1991; 75 Lawson et al. 2008; Saeed et al. 2022; Ud Din et al. 2018).

International Standards Organisation (ISO) 4823 (2007), recommends that an elastomeric impression material (all consistencies) must have 96.5% elastic recovery. Although all elastomeric impression materials fulfil the criteria, the VPS possesses better elastic recovery compared to other impression materials (Bonsor & Pearson 2013; Din et al. 2018b; Hamalian et al. 2011; Klooster et al. 1991). This allows pouring of the impression to fabricate cast after six minutes of removal from the mouth (Blomberg et al. 1992). In addition, these materials exhibit great dimensional stability and high tear strength.

84

Different brands of VPS impression materials showed variations in elastic recovery. Lawson et al. (2008) investigated the elastic recovery for five VPS and a hybrid impression material (containing siloxane and polyether groups) after subjecting materials to tensile and compressive stress. The VPS impression materials showed improved elastic recovery in comparison to the hybrid material, which may be related to the compositions of materials as hybrid material composed of polyethers, which have a lower elastic recovery compared to VPS (Hondrum 1994; Ud Din et al. 2022a; Ud Din et al. 2022b).

92

93 Strain-in-compression is the measurement of the stiffness or flexibility of impression material. It 94 determines the ability of polymerized material to be removed from the mouth or cast without 95 permanent deformation, injury to oral tissues or fracture. Also, it dictates the ability of the 96 impression to resist deformation and withstand the weight of the dental stone when the cast is 97 poured (Helvey 2011; Klooster et al. 1991; Lu et al. 2004a).

99 To overcome the problem of inherent hydrophobicity of VPS and to improve tear strength and % 100 elongation at break of the material, in our previous work, novel formulations of VPS were 101 fabricated using a novel cross-linking agent i.e. tetra-functional (dimethylsilyl) ortho-silicate 102 (TFDMSOS) and novel surfactant i.e. Rhodasurf CET-2 (ethoxylated cetyloleyl alcohol (Din et al. 103 2018a; Din et al. 2018b; Din et al. 2017). The addition of TFDMOS improved mechanical 104 properties of experimental impression materials, while the non-ionic surfactant was added to 105 improve wetting properties of the materials. Different researchers have explored the effect of 106 various surfactants to improve the hydrophilicity of the material, however, little work has been 107 carried out to improve the tear strength of VPS impressions. Additionally, the effects of the 108 addition of surfactant on the mechanical properties of the materials and the hydrophilicity of these 109 modified materials after disinfection requires further exploration. Rhodasurf CET-2 is a non-ionic 110 surfactant which is made by combination of ethoxylated cetyl and ethoxylated oleyl alcohols. 111 Ethoxylated oleyl alcohol has a double bond in its chemical structure. The double bond could be 112 possibly activated during mixing and played a vital role in the cross-linking polymerization 113 reaction leading to improved elastic recovery and strain-in-compression²³.

114 Ud Din et al. observed that the incorporation of a novel cross-linking agent (TFDMSOS) 115 significantly improved the materials' % elongation-at-break and tear strength compared to the 116 control containing a conventional crosslinking agent-poly (methyl-hydro siloxane)⁵. Additionally, 117 the incorporation of a novel surfactant (Rhodasurf CET-2) further resulted in a significant increase 118 in % elongation-at-break (Din et al. 2018a). It was also noted that the experimental formulation 119 has a lower contact angle (improved hydrophilicity) than commercial formulations. Additionally, 120 disinfection has little impact on the contact angle as the surfactant did not readily leach out in a 121 disinfecting solution (Din et al. 2017). However, mechanical testing including elastic recovery and strain-in-compression required further exploration before considering the experimentalformulation as a better alternative to commercially available VPS impression materials.

124

125 The purpose of this study was to evaluate the effect of a novel cross-linking agent, TFDMSOS and 126 novel surfactant (Rhodasurf CET-2) on the elastic recovery and strain-in-compression of 127 experimental VPS and to compare it with commercial materials. In summary, addition silicone 128 materials with higher cross-link density have better elastic recovery, as they have a greater number 129 of cross-links that can resist deformation. This property, along with other desirable characteristics 130 such as tear strength and dimensional stability, make addition silicone impression materials a 131 popular choice for dental impressions. It was hypothesized that in the current study the 132 experimental formulations have better elastic recovery and strain-in compression-values due to 133 higher cross-links provided by adding TFDMSOS and Rhodasurf CET-2 and making it a more 134 suitable material for recording an accurate impression.

135

136 Materials & Methods

The ethical approval and informed consent were not required for this study, since this study do not involve living human subjects and only involve *in vitro* laboratory testing of dental impression material. Three medium-body commercial VPS impression materials were used in this study; Aquasil Ultra Monophase, USA, Dentsply (Aq M), Elite HD Monophase, Italy, Zhermack (Elt M) and Extrude, USA, Kerr (Extr M). Additionally, five experimental VPS formulations were prepared as base paste and catalyst paste (Table 1) ⁵. Exp-I was used as a control for Exp-II, while Exp-II acted as a control group for Exp-III, Exp-IV and Exp-V.

145 **Preparation of Experimental Formulations**

The base paste of Exp-1 (hydrophobic VPS) was formulated by mixing vinyl-terminated polydimethyl siloxane and a conventional cross-linking agent (poly methylhydrosiloxane) for 5 minutes using an electric hand mixer (Kenwood, kMix, UK). The filler (Aerosil R812S) was added to the mixture and a uniform paste was made by mixing the components with a pestle and mortar for 5 min, followed by blending the paste with an electric mixer for 10 minutes.

151

The catalyst paste was formulated by mixing vinyl-terminated poly (dimethylsiloxane), platinum catalyst and palladium for 5 minutes with the electric hand mixer, followed by the addition of filler (Aerosil R812S) and mixing it with the pestle and mortar and electric hand mixer. For preparation of Exp-II impression material, the amount of poly (methylhydrosiloxane) was reduced from 1.10% to 0.77% and it was replaced it with a novel cross-linking agent (TFDMSOS) in the base paste of Exp-I formulation. Vinyl to silane groups were maintained at 1:1 ratio. The catalyst paste for Exp-II was similar to that of Exp-I.

Experimental formulations III, IV and V were formulated by modifying base-paste of Exp-II with addition of non-ionic surfactant; Rhodasurf CET-2) at concentrations of 2%, 2.5% and 3% respectively. The quantities of constituents in the catalyst paste were adjusted to ensure adequate polymerization of the materials (Table 1). The prepared base and catalyst paste of experimental materials were kept in separate compartments of cartridge and stored at 4°C \pm 2°C.

164 Sample preparation for elastic recovery and strain under compression

The cylindrical samples for elastic recovery and strain-in-compression were prepared using polytetrafluoroethylene (PTFE) mould measuring 20 mm in length x 12.5 mm in diameter according to ISO4823 (2007) standard. PTFE mould was positioned on top of a metal plate lined with an acetate sheet. The base and catalyst pastes were syringed into the mould using an automixing syringe and the mould was sandwiched by another acetate-lined metal plate. The assembly was held using C-clamp. Commercial materials were left to polymerize according to the manufacturer's instructions while experimental materials were allowed to be set for 4 to 11 minutes (Din et al. 2017).

173 To measure elastic recovery (n=12), two metal plates (13 x 13 x 3 mm³) were fixed on either side 174 of the specimen with the aid of double-sided sticky tape. The length of the specimen including 175 metal plates (h₁) was recorded using a digital micrometre (Mitutoyo, Japan) to an accuracy of 176 0.001 mm. Then the specimen was deformed to 6 ± 0.1 mm within 1 second using the calibrated 177 mechanical testing machine (Tinius Olsen Ltd, Model H5KS, England, load cell 5kN) shown in 178 Figure 1. The deformation force was released slowly over a period of 5 seconds. After two 179 minutes the length was measured again (h₂). The elastic recovery in percentage, K, was assessed 180 using equation 1.

181

182 Equation 1: K = 100 -
$$\left[100\left(\frac{h_1 - h_2}{ho}\right)\right]$$

183 h_0 is the height (mm) of the split mould

184 h₁ is the length (mm) of the specimen immediately before the application of the initial load

185 h₂ is the length of the specimen, 2 minutes after removing the deformation force

186

To evaluate strain-in-compression, 12 samples per material were tested. An initial force of
1.22±0.1 N was exerted on the specimen and the distance (h₁) was calculated using the Tinius
Olsen (Figure 1). The load was increased to 12.25±0.1 (N) progressively over a time of 10 seconds

at a rate of 3 mm/min and a change in height of the specimen was noted again (h₂). The percentage
of strain-in-compression, E, was calculated using equation 2.

192

193 Equation 2:
$$E = \left(\frac{h_1 - h_2}{h_0}\right) 100$$

194 ho is the height (mm) of the split mould

195 h₁ is the length (mm) of the specimen, 30 seconds after submission of the opening load

 h_2 is the length of the specimen, 30 seconds after submission of the amplified load.

197

198 The data was analyzed using SPSS Version 22 (Armonk NY IBM Corp, Armonk, NY, USA).

199 Numerical data were presented as mean and standard deviation. Analysis of variance was

200 performed with p value at 0.05. Where significant difference in group was found, individual

201 means were compared using post hoc Tukey's test.

202

203 **Results**

Table 2 shows the elastic recovery and strain-in-compression for commercial and experimental

205 VPS impression materials immediately after setting. All the tested materials met the ISO4823

requirement of having elastic recovery greater than 96.5%. Exp-II exhibited the highest elastic

207 recovery while Exp-V demonstrated the lowest values. The post-hoc analysis revealed that all

208 three commercial products and Exp-I had statistically similar elastic recoveries.

209

Use of novel crosslinking agent (TFDMSOS) instead of conventional agent significantly increased
elastic recovery. A significant difference in the elastic recovery was noted between Exp-II and
Exp-V. It was noted that the addition of a non-ionic surfactant (Rhodasurf CET-2) in the

experimental formulation, to improve hydrophilicity of material, resulted in a reduced percentageof elastic recovery of material, however, the changes were statistically not significant (Table 2).

215

216 Strain-in-compression

217 Figure 2 and Table 2 reveal the strain-in-compression for the tested VPS impression materials. 218 Experimental VPS had significantly higher (p <0.05) strain-in-compressions values compared to 219 the commercial VPS. Exp-V exhibited significantly the highest (Tukey's HSD test) strain-in-220 compression (7.08 % \pm 0.22 %) while Elt M had the lowest values (3.15 % \pm 0.18 %). Among 221 commercial materials, no significant difference between Aq M and Extr M was noted. However, it was noted that the addition of a novel crosslinking agent i.e., TFDMSOS (Exp II), significantly 222 223 increased the percentage strain-in-compression values compared to formulations based on 224 conventional cross-linking agents (Exp-I, Aq M, Elt M, Extr M). Also, it was observed that 225 experimental formulations incorporating non-ionic surfactant (Rhodasurf CET-2) led to a further 226 significant increase in strain-in-compression values and this effect was concentration dependant.

227

228 Figure 3 and Table 2 show the comparison between elastic recovery and strain-in-compression for 229 all commercial and Experimental VPS impression materials evaluated in this study. Among the 230 experimental materials, there is a correlation between elastic recovery and strain-in-compression. 231 With the addition of TFDMSOS in Exp-II the elastic recovery and strain-in-compression increase 232 significantly compared to Exp-I (control). However, there is a negative correlation seen after 233 addition of Rhodasurf CET-2 (non-ionic surfactant) in Exp III. With the addition of surfactant the 234 elastic recovery is decreased while stain-in-compression is increased. It can also be noticed that 235 with the increase in % amount of surfactant there is a consistent and significant decrease in elastic

recovery and significant incurease in strain-in-compression in Exp-IV and Exp-V. Amongcommercial materials, no significant difference was seen.

238

239 **Discussion**

The elastic recovery of the impression materials plays a major role in the accurate reproduction of the oral cavity. The ability of elastomeric impression materials to revert to their actual form upon removal of the applied stress is related to their coiled wrapped and kinked molecular chains. These polymer chains stretch in the direction of stress and, recoil back on releasing the force, gaining their original shape and form (Hamalian et al. 2011; Klooster et al. 1991). In the present study, elastic recovery and strain-in-compression of commercial and experimental VPS impression materials were compared.

247 Values for elastic recovery for commercial and experimental VPS impression materials ranged 248 from 99.11 to 99.64%. These values were within the range set by International Standards 249 Organization (ISO) 4823 (2007) which requires \geq 96.5 %. Similar results were reported by Lawson 250 et al. (2008), who noted that elastic recovery of five tested VPS (Aquasil Ultra, Examix, Genie, 251 Imprint 3, and StandOut) and one hybrid impression material (Senn) ranges from 99.34 to 99.83 252 %. In another study Lu et al. (2004a) investigated the elastic recovery of two VPS (Flexitime and 253 Imprint II) and one polyether (Impregum). It was noted that Flexitime, Imprint II and Impregum 254 had 99.60, 99.75 and 99.19 % elastic recoveries respectively.

Similar percentage elastic recovery was noted for the commercial materials (Aq M: 99.32%, Elt
M: 99.31%, Extr M: 99.27%) and Exp-I (99.32%) containing conventional cross-linking agent
polymethylhydrosiloxane. However, on incorporating a novel cross-linking agent (TFDMSOS) in
Exp-II, an increase in elastic recovery (99.65%) was observed which was statistically significant.

259 The greater elastic recovery of Exp-II is attributed to excellent crosslinking of TFDMSOS with 260 functional groups of vinyl-terminated poly (dimethylsiloxane) pre-polymer as each molecule of TFDMSOS can bond to four functional groups of pre-polymer (Din et al. 2018a). Similar results 261 262 have been reported in the literature indicating the amount of permeant deformation of an 263 impression material is greatly influenced by the degree of cross-linking of the polymeric chains 264 (Singer et al. 2022). The degree of polymerization also affects other mechanical properties of 265 elastomeric impression materials such as tear strength and % elongation-at-break (Din et al. 2018a; 266 Din et al. 2018b). It is reported that both the properties improved due to increased cross-linking of 267 polymeric chains.

268 The strain-in-compression was also calculated to assess the rigidity of impression materials so that 269 it can be removed from the mouth or cast without permanent deformation after setting, and to resist 270 deformation when the dental stone is poured. All tested impression materials have values for % 271 strain-in-compression within the ISO4823 (2007) limits. Experimental VPS impression materials 272 had higher strain-in-compression values indicating improved flexibility of experimental material 273 (Lu et al. 2004b). Therefore, a positive correlation between elastic recovery and strain-in-274 compression was noted (Figure 3 and Table 2). Additionally, it was observed that the incorporation 275 of the wetting agent (Rhodasurf CET-2), further significantly increased strain-in-compression 276 values. This was contradictory to the results of Lu et al. who noticed that flexible materials have 277 less cross-linking and have better elastic recovery. This might be due to the difference in the 278 composition of the materials used in the present study.

The ability to undergo greater elastic recovery is a desirable property of impression materials as it ensures an accurate impression which in turn guarantees a correct fit of the prosthesis. The experimental VPS impression materials in this study show greater elastic recovery than their

283	proven to have improved wettability, percentage elongation, tear strength and minimal distortion			
284	after disinfection making them a much more suitable option for impression taking (Din et al.			
285	2018a; Din et al. 2018b; Din et al. 2021).			
286	A limitation of this study is that it is conducted in an in-vitro environment under laboratory			
287	conditions. To strengthen the claim of experimental VPS as a superior impression material, it is			
288	necessary to conduct further research in intra-oral, in-vivo, conditions.			
289				
290	Conclusions			
291	The addition of a novel cross-linking agent (TFDMSOS) showed improved elastic recovery and			
292	strain-in-compression, while the addition of a non-ionic surfactant also significantly increased			
293	strain-in-compressions values for all experimental VPS. All tested materials comply with ISO			
294	standards. In the future, biocompatibility testing followed by clinical trials should be undertaken,			
295	and material selection should be based on adequate knowledge of the properties of materials to			
296	improve clinical success.			
297				
298				
299	References			
300	Balkenhol M, Haunschild S, Erbe C, and Wöstmann B. 2010. Influence of prolonged setting time on			
301	permanent deformation of elastomeric impression materials. The Journal of Prosthetic Dentistry			
302	103:288-294.			
303	Blomberg PA, Mahmood S, Smales RJ, and Makinson OF. 1992. Comparative elasticity tests for			
304	elastomeric (non putty) impression materials. <i>Australian Dental Journal</i> 37:346-352.			
305	Bonsor SJ, and Pearson GJ. 2013. A Clinical Guide to Applied Dental Materials Churchill Livingstone			
300	Elseviel. Braden M Clarke RI Nicholson I and Parker S 1997 Polymeric Dental Materials Springer- Verlag			
308	Darvell WB 2009 Materials Science for Dentistry 9th ed: B W Darvell			
309	Deb S. 1998. Polymers in dentistry. <i>Proceedings of the Institution of Mechanical Engineers</i> . Part H:			
310	Journal of Engineering in Medicine 212:453-464.			
311	Din SU, Noor N, Humayoun S, Khalid S, Parker S, and Patel M. 2018a. Tensile strength of novel			
312	experimental hydrophilic vinyl polysiloxane impression materials compared to control and			

commercial counterparts. Additionally, in previous studies, same experimental material has

282

313 commercial VPS impression materials. *Journal of Islamabad Medical & Dental College* 7:67-72.

- Din SU, Parker S, Braden M, and Patel M. 2018b. The effects of cross-linking agent and surfactant on the
 tear strength of novel vinyl polysiloxane impression materials. *Dental materials* 34:e334-e343.
- Din SU, Parker S, Braden M, and Patel M. 2021. Improved water absorption behaviour of experimental
 hydrophilic vinyl polysiloxane (VPS) impression materials incorporating a crosslinking agent and
 a novel surfactant. *Dental materials* 37:1054-1065.
- Din SU, Parker S, Braden M, Tomlins P, and Patel M. 2017. Experimental hydrophilic vinyl polysiloxane
 (VPS) impression materials incorporating a novel surfactant compared with commercial VPS.
 Dental materials 33:e301-e309.
- Din SU, Sajid M, Saeed A, Chaudhary FA, Alam MK, Sarfraz J, Ahmed B, and Patel M. 2022.
 Dimensional changes of commercial and novel polyvinyl siloxane impression materials following sodium hypochlorite disinfection. *PeerJ* 10:e12812.
- Goldberg AJ. 1974. Viscoelastic properties of silicone, polysulphide, and polyether impression materials.
 Journal of Dental Research 53:1033-1039.
- Hamalian TA, Nasr E, and Chidiac JJ. 2011. Impression materials in fixed prosthodontics: Influence of
 choice on clinical procedure. *Journal of Prosthodontics* 20:153-160.
- Helvey GA. 2011. Elastomeric impression materials: factors to consider. *Compendium of continuing education in dentistry (Jamesburg, NJ : 1995)* 32:58-59.
- Hondrum SO. 1994. Tear and energy properties of three impression materials. *The International journal of prosthodontics* 7:517-521.
- 333 ISO4823. 2007. Dentistry Elastomeric Impression Materials ISO 4823:2000/Amd.1.2007(E). p 1-33.
- Klooster J, Logan GI, and Tjan AHL. 1991. Effects of strain rate on the behavior of elastomeric
 impression. *The Journal of Prosthetic Dentistry* 66:292-298.
- Lawson NC, Burgess JO, and Litaker MS. 2008. Tensile elastic recovery of elastomeric impression
 materials. *Journal of Prosthetic Dentistry* 100:29-33.
- Lu H, Nguyen B, and Powers J. 2004a. Mechanical properties of 3 hydrophilic addition silicone and
 polyether elastomeric impression materials* 1. *The Journal of Prosthetic Dentistry* 92:151-154.
- Lu H, Nguyen B, and Powers JM. 2004b. Mechanical properties of 3 hydrophilic addition silicone and
 polyether elastomeric impression materials. *The Journal of prosthetic dentistry* 92:151-154.
- 342 Mandikos M. 1998. Polyvinyl siloxane impression materials: An update on clinical use. *Australian* 343 *Dental Journal* 43:428-434.
- Saeed A, Zahid S, Sajid M, Ud Din S, Alam MK, Chaudhary FA, Kaleem M, Alswairki HJ, and
 Abutayyem H. 2022. Physico-Mechanical Properties of Commercially Available Tissue
 Conditioner Modified with Synthesized Chitosan Oligosaccharide. *Polymers* 14:1233.
- Singer L, Bourauel C, Habib SI, Shalaby HE-A, and Saniour SH. 2022. Tear strength and elastic recovery
 of new generation hybrid elastomeric impression material: A comparative study. *BMC Research Notes* 15:1-5.
- Ud Din S, Chaudhary FA, Ahmed B, Alam MK, Parker S, Patel M, and Javed MQ. 2022a. Comparison of
 the Hardness of Novel Experimental Vinyl Poly Siloxane (VPS) Impression Materials with
 Commercially Available Ones. *BioMed research international* 2022.
- Ud Din S, Chaudhary FA, Alyahya Y, Arjumand B, Javed MQ, Ali A, and Patel M. 2022b. Reproduction
 of Fine Details and Compatibility of Vinyl Polysiloxane Impression Materials. *Coatings* 12:867.
- Ud Din S, Hassan M, Khalid S, Zafar MS, Ahmed B, and Patel M. 2018. Effect of surfactant's molecular
 weight on the wettability of vinyl polysiloxane impression materials after immersion disinfection.
 Materials Express 8:85-92.
- 358





Figure 1: Tinius Olsen





Commente		Base Paste (Wt %)				
Components	Exp-I	Exp-II	Exp-III	Exp-IV	Exp-V	
Vinyl-terminated dimethylpolysiloxane (Mw 62700)	39.90	39.90	37.95	37.46	36.98	
Polymethylhydrosiloxane (Mw 2270)	1.10	0.77	0.74	0.73	0.72	
Tetra-functional (dimethylsilyl) orthosilicate (TFDMSOS) (Mw 329)	-	0.33	0.32	0.31	0.31	
Filler Aerosil R 812	9	9	9	9	9	
Components		Catalyst Paste (Wt %)				
	Exp-I	Exp-II	Exp-III	Exp-IV	Exp-V	
Vinyl-terminated dimethylpolysiloxane (Mw 62700)	40.72	40.72	39.51	39.51	39.51	
Platinum (0.05 M)	0.06	0.06	1.27	1.27	1.27	
Palladium (<1µm)	0.23	0.23	0.22	0.22	0.22	
Filler Aerosil R 812	9	9	9	9	9	
Rhodasurf CET-2 (non-ionic surfactant)	-	-	2.00	2.50	3.00	

Table 1. Composition of novel experimental (Exp-I, II, III, IV and V) VPS impression materials.

Table 2. Average elastic recovery and strain in compression of commercial and experimental VPS immediately after setting.

Impression Materials	Elastic Recovery (%)	Strain-in-compression
Aq M	99.32 ± 0.30	4.261±0.154
Elt M	99.31 ± 0.35	3.153±0.177
Extr M	99.27 ± 0.32	4.405±0.118
Exp-I	99.32 ± 0.52	4.677±0.207
Exp-II	99.65 ± 0.09	5.360±0.163
Exp-III	99.50 ± 0.23	6.137±0.256
Exp-IV	99.37 ± 0.26	6.541±0.239
Exp-V	99.12 ± 0.16	7.076±0.220